

# Motorized Examination Solutions

## Mi-4004 B



\*The above product can be supplied with either a foot switch or a hand control

### Specifications

Length x Width	1875 mm x 615 mm
Minimum/Maximum Height	500 mm (20") / 800 mm (32")
Patient Weight Bearing Capacity	150 kg
Backrest Angle Adjustment	0-65°

### Standard Product Features

Backrest and Height Functions electrically operated  
Paper Roll Holder  
Single Drawer at the Foot End of the Table

### Optional Features

Battery Back up  
Choice of Hand Control/ Foot Switch